



# SSP-02-ZENIT


MEDICAL PORTABLE LAMP

SCU#: 36703

SSP-02-Zenit medical portable lamp is developed for lighting the working field during surgical procedures and examinations in clinics and hospitals

## SSP-02-ZENIT

With LED lamps and fixation unit for attachment to a medical chair



### SPECIFICATIONS:

Working distance, m	0.8
Working field diameter, mm	
• LED lamps	> 200
• halogen lamps	—
Maximum luminance in the center of the working field at a distance of 0.8 m, lux	
• LED lamps ( $\pm 10000$ lc)	40000
• halogen lamps ( $\pm 2500$ lc)	—
Specific irradiance of the working field, $W\ m^{-2}/klk$	< 5.0
Power consumption, V*A	< 35
Power	AC single-phase network with frequency $50 \pm 0.5$ Hz, Voltage, V: $220 \pm 22$
Dimensions (HxWxD),mm	1000 x 200 x 150
Weight, kg	$1.8 \pm 0.5$



Rostec



Shvabe



PJSC "Krasnogorsky Zavod"

Tel.: +7 (495) 561-27-08, +7 (495) 561-80-08

E-mail: [gp@zenit-kmz.ru](mailto:gp@zenit-kmz.ru)

[www.zenit-medicine.com](http://www.zenit-medicine.com)

[www.shvabe.com](http://www.shvabe.com)