



SSP-02-6/1 ZENIT

MEDICAL PORTABLE LAMP

SSP-02-6/1 Zenit medical portable lamp is developed for lighting the working field during surgical procedures and examinations in clinics and hospitals

# SSP-02-6/1 ZENIT

With the ring LED lighting unit, the central halogen lamp and the five-wheel base support without brake



## SPECIFICATIONS:

Working distance, m	0.8
Working field diameter, mm	
<ul style="list-style-type: none"> <li>• LED lamps &gt; 200</li> <li>• halogen lamps &gt; 100</li> </ul>	
Maximum luminance in the center of the working field at a distance of 0.8 m, lux	
<ul style="list-style-type: none"> <li>• LED lamps (<math>\pm 10000</math> lc) 40000</li> <li>• halogen lamps (<math>\pm 2500</math> lc) 7500</li> </ul>	
Specific irradiance of the working field, $W\ m^{-2}/klk$	< 5.0
Power consumption, V*A	< 95
Power	AC single-phase network with frequency $50 \pm 0.5$ Hz, Voltage, V: $220 \pm 22$
Dimensions (HxWxD),mm	1800 x 600 x 600
Weight, kg	$5 \pm 0.5$



Rostec



Shvabe



PJSC "Krasnogorsky Zavod"

Tel.: +7 (495) 561-27-08, +7 (495) 561-80-08

E-mail: [gp@zenit-kmz.ru](mailto:gp@zenit-kmz.ru)

[www.zenit-medicine.com](http://www.zenit-medicine.com)

[www.shvabe.com](http://www.shvabe.com)