




# LS-01 – ZENIT

## SLIT LAMP

- W/O DIGITAL CAMERA
- EQUIPPED WITH OPHTHALMIC TABLE

SCU#: 39001



## LS-01–ZENIT

Excellent optical characteristics of the lamp are obtained due to the use of Grenou optical system, which provides a real stereo effect, a greater depth of field, high resolution, a better ratio of zooming to the covered field and a convergent beam path

### SPECIFICATIONS:

Type of stereo microscope	Grenou
Three level zooming, x	8, 16, 32
Working field diameter, mm	28 , 14.5, 7.5
Ocular adjustment, dptr	-7 to +7
Interpupillary distance adjustment, mm	58 to 75
Filters	Blue, Green, Neutral, IR



Rostec



Shvabe



ZENIT

PJSC "Krasnogorsky Zavod"

Tel.: +7 (495) 561-27-08, +7 (495) 561-80-08

E-mail: [gp@zenit-kmz.ru](mailto:gp@zenit-kmz.ru)

[www.zenit-medicine.com](http://www.zenit-medicine.com)

[www.shvabe.com](http://www.shvabe.com)

### DELIVERY SET:

- Slit lamp LS-01–«Zenit»
- Ophthalmic table with electric drive