



KM-01-ZENIT

GYNECOLOGICAL CHAIR

SCU#: 371

KM-01-ZENIT

Zenit KM-01 gynecological chair with three electric motors is the most reliable model among competitive products. Equipped with independent electromechanical motors the chair provides the convenient positioning of a patient during medical examinations

DESIGN AND ERGONOMICS FOR INTELLIGENT WORKING CONVENIENCE:

- two built-in sockets for connecting medical lamps and other equipment;
- Trendelenburg/Anti-Trendelenburg position adjustment of the seat bed;
- 3 motors for height, backrest and Trendelenburg position adjustment;
- the installed universal guide rails provide a convenient attachment of accessories on the chair;
- Goepel-type leg supports with attachment clamps;
- a removable stainless steel bowl

DURABLE PERFORMANCE DUE TO:

- LINAK electric motors;
- Arpatek eco-leather upholstery that withstands 500,000 cycles according to the Martindale test, that is 16 times more than required for heavy domestic conditions or frequent commercial use

TECHNICAL SPECIFICATIONS:

Bed width, mm	600
Bed length, mm	1250
Minimum seat height, mm	880
Load capacity, kg	140
Total weight (with accessories), kg	<180

MOTION PARAMETERS (THREE ELECTRIC MOTORS):

Overall height adjustment, mm	880 – 1085
Seat position adjustment (to the backrest)	135° - 180°
Tilt angle adjustment between back and seat	
Bed position adjustment (back + seat)	- 10° (down) - + 45° (up)



Rostec



Shvabe



PJSC "Krasnogorsky Zavod"

Tel.: +7 (495) 561-27-08, +7 (495) 561-80-08

E-mail: gp@zenit-kmz.ru

www.zenit-medicine.com

www.shvabe.com