



KNB-02-ZENIT

BINOCULAR FLOOR COLPOSCOPE

- HALOGEN LAMP

SCU#: 13303

KNB-02-ZENIT BINOCULAR FLOOR COLPOSCOPE

KNB-02 -Zenit colposcope is designed for medical monitoring during examinations and diagnostics of gynecological diseases and in obstetrics via zooming a stereoscopic image of the working field

KEY POINTS :

- Working distance: 300 mm (optimal);
- Converged optical system reduces eye fatigue during routine examinations;
- 3D image allows to view the structure and condition of blood vessels in detail and the blue/green filter helps to select them in contrast;
- The integrated tripod provides rotation of the stereo microscope around the axis at any angle, adjustment of height and focus.



SPECIFICATIONS:

Working distance, mm	300
Zooming, x	7.5; 15; 30
Working field diameter, mm	50; 25; 12
Interpupillary distance adjustment, mm	58~75
Dioptric adjustment, dptr	-7 to +7

Light system:

Light source	halogen lamp
Maximum luminance (adjustable), lux	> 10000
Covered field diameter, mm	> 70



Shvabe



PJSC "Krasnogorsky Zavod"

Tel.: +7 (495) 561-27-08, +7 (495) 561-80-08

E-mail: gp@zenit-kmz.ru

www.zenit-medicine.com

www.shvabe.com